

**NUCLEODUR C18 Isis, 3 µm****Macherey Nagel****NUCLEODUR® C18 Isis, 3 µm**

Fase C18 con polimeri speciali, USP L1

dimensioni pori 110Å, dimensioni particelle 1,8 µm, 3 µm e 5 µm; 20% C

alta selettività sterica

Notevole superficie di deattivazione

idonea per LC/MS grazie alla bassa capacità legante.

stabilità pH da 1 a 10

ampio range di applicazioni: steroidi, (o,p,m), vitamine solubili ect.

colonne analitiche EC, dimensioni particelle 3 µm, 20% C

Ø int. 2.0 mm

Lungh. 50 mm



Code	Description	Packaging
MNL76040020	NUCLEODUR C18 Isis, 3 µm, Ø int. 2.0 mm, Lungh. 50 mm	1 pz
MNL76040030	EC 50/3 NUCLEODUR C18 Isis, 3 µm EC HPLC column EC 50/3 NUCLEODUR C18 Isis, 3 µm length: 50 mm, ID: 3 mm	1 pz
MNL76040040	NUCLEODUR C18 Isis, 3 µm, Ø int. 4.0 mm, Lungh. 50 mm	1 pz
MNL76040046	EC 50/4.6 NUCLEODUR C18 Isis, 3 µm EC HPLC column EC 50/4.6 NUCLEODUR C18 Isis, 3 µm length: 50 mm, ID: 4.6 mm	1 pz
MNL76040120	EC 100/2 NUCLEODUR C18 Isis, 3 µm EC HPLC column EC 100/2 NUCLEODUR C18 Isis, 3 µm length: 100 mm, ID: 2 mm	1 pz
MNL76040130	EC 100/3 NUCLEODUR C18 Isis, 3 µm EC HPLC column EC 100/3 NUCLEODUR C18 Isis, 3 µm length: 100 mm, ID: 3 mm	1 pz
MNL76040140	EC 100/4 NUCLEODUR C18 Isis, 3 µm EC HPLC column EC 100/4 NUCLEODUR C18 Isis, 3 µm length: 100 mm, ID: 4 mm	1 pz
MNL76040146	NUCLEODUR C18 Isis, 3 µm, Ø int. 4.6 mm, Lungh. 100 mm	1 pz
MNL76040220	NUCLEODUR C18 Isis, 3 µm, Ø int. 2.0 mm, Lungh. 125 mm	1 pz
MNL76040230	NUCLEODUR C18 Isis, 3 µm, Ø int. 3.0 mm, Lungh. 125 mm	1 pz
MNL76040240	NUCLEODUR C18 Isis, 3 µm, Ø int. 4.0 mm, Lungh. 125 mm	1 pz
MNL76040246	NUCLEODUR C18 Isis, 3 µm, Ø int. 4.6 mm, Lungh. 125 mm	1 pz
MNL76040320	NUCLEODUR C18 Isis, 3 µm, Ø int. 2.0 mm, Lungh. 150 mm	1 pz
MNL76040330	NUCLEODUR C18 Isis, 3 µm, Ø int. 3.0 mm, Lungh. 150 mm	1 pz
MNL76040340	NUCLEODUR C18 Isis, 3 µm, Ø int. 4.0 mm, Lungh. 150 mm	1 pz
MNL76040346	NUCLEODUR C18 Isis, 3 µm, Ø int. 4.6 mm, Lungh. 150 mm	1 pz
MNL76040420	NUCLEODUR C18 Isis, 3 µm, Ø int. 2.0 mm, Lungh. 250 mm	1 pz
MNL76040430	NUCLEODUR C18 Isis, 3 µm, Ø int. 3.0 mm, Lungh. 250 mm	1 pz



Code	Description	Packaging
MNL76040440	NUCLEODUR C18 Isis, 3 µm, Ø int. 4.0 mm, Lungh. 250 mm	1 pz
MNL76040446	NUCLEODUR C18 Isis, 3 µm, Ø int. 4.6 mm, Lungh. 250 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

