

**NUCLEODUR Sphinx RP, 5 µm****Macherey Nagel****NUCLEODUR® Sphinx RP, 5 µm**

selettività distintiva basata sulla superficie bifunzionale - USP L1 e USP L11
dimensione pori 110Å, dimensione particelle 1,8 µm, 3 µm e 5 µm; 14 % C
alta densità di legami covalenti silani per picchi scodati per allargare lo sviluppo del metodo
stabilità pH 1 - 10
adatte per LC/MS grazie alle caratteristiche di basso spurgo
alta riproducibilità e consistenete qualità grazie alle strette procedure QC
range applicativo: antibiotici al chinolone, sulfamidici, xantine, sostituti aromatici
colonne analitiche EC, dimensione particelle 5 µm, 14% C

Ø int. 2.0 mm
Lungh. 50 mm



Code	Description	Packaging
MNL76080020	NUCLEODUR Sphinx RP, 5 µm, Ø int. 2.0 mm, Lungh. 50 mm	1 pz
MNL76080030	EC 50/3 NUCLEODUR SPHINX RP, 5 µm EC HPLC column EC 50/3 NUCLEODUR SPHINX RP, 5 µm length: 50 mm, ID: 3 mm	1 pz
MNL76080040	EC 50/4 NUCLEODUR SPHINX RP, 5 µm EC HPLC column EC 50/4 NUCLEODUR SPHINX RP, 5 µm length: 50 mm, ID: 4 mm	1 pz
MNL76080046	EC 50/4.6 NUCLEODUR SPHINX RP, 5 µm EC HPLC column EC 50/4.6 NUCLEODUR SPHINX RP, 5 µm length: 50 mm, ID: 4,6 mm	1 pz
MNL76080120	EC 125/2 NUCLEODUR SPHINX RP, 5 µm EC HPLC column EC 125/2 NUCLEODUR SPHINX RP, 5 µm length: 125 mm, ID: 2 mm	1 pz
MNL76080130	EC 125/3 NUCLEODUR SPHINX RP, 5 µm EC HPLC column EC 125/3 NUCLEODUR SPHINX RP, 5 µm length: 125 mm, ID: 3 mm	1 pz
MNL76080140	NUCLEODUR Sphinx RP, 5 µm, Ø int. 4.0 mm, Lungh. 125 mm	1 pz
MNL76080146	NUCLEODUR Sphinx RP, 5 µm, Ø int. 4.6 mm, Lungh. 125 mm	1 pz
MNL76080220	EC 150/2 NUCLEODUR SPHINX RP, 5 µm EC HPLC column EC 150/2 NUCLEODUR SPHINX RP, 5 µm length: 150 mm, ID: 2 mm	1 pz
MNL76080230	EC 150/3 NUCLEODUR SPHINX RP, 5 µm EC HPLC column EC 150/3 NUCLEODUR SPHINX RP, 5 µm length: 150 mm, ID: 3 mm	1 pz
MNL76080240	NUCLEODUR Sphinx RP, 5 µm, Ø int. 4.0 mm, Lungh. 150 mm	1 pz
MNL76080246	EC 150/4.6 NUCLEODUR SPHINX RP, 5 µm EC HPLC column EC 150/4.6 NUCLEODUR SPHINX RP, 5 µm length: 150 mm, ID: 4,6 mm	1 pz
MNL76080320	NUCLEODUR Sphinx RP, 5 µm, Ø int. 2.0 mm, Lungh. 250 mm	1 pz
MNL76080330	EC 250/3 NUCLEODUR SPHINX RP, 5 µm EC HPLC column EC 250/3 NUCLEODUR SPHINX RP, 5 µm length: 250 mm, ID: 3 mm	1 pz
MNL76080340	NUCLEODUR Sphinx RP, 5 µm, Ø int. 4.0 mm, Lungh. 250 mm	1 pz





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

