


**VarioPrep supporti per precolonne e parti di ricambio****Macherey Nagel****VarioPrep supporti per precolonne e parti di ricambio
colonne VarioPrep**

sistema colonna per HPLC preparativa prodotto in acciaio inossidabile con un terminale regolabile (su richiesta sono disponibili anche colonne con due terminali regolabili, ad esempio per uso frequente di tecniche di riflusso) consente compensazione del volume morto, che potrebbe presentarsi all'ingresso della colonna dopo diversi utilizzi, senza necessità di aprire la colonna impaccata con silice sferica NUCLEODUR® e NUCLEOSIL®. O-rings disponibili a richiesta

Descrizione Supporto VP per precolonne 8mm per colonne VarioPrep con diam.int. 8 e 10mm



Code	Description	Packaging
MNL718251	VarioPrep supporti per pecolonne e parti di ricambio, Supporto VP per precolonne 8mm per colonne VarioPrep con diam.int. 8 e 10mm	1 pz
MNL718253	VP Guard column holder 32 mm (short) VarioPrep guard column holder (short) for 32 mm ID guard column cartridges VP 15/32 recommended for VarioPrep columns with 32 and 40 mm ID	1 pz
MNL718255	VP guard column holder 50 mm (short) VarioPrep guard column holder (short) for 50 mm ID guard column cartridges VP 15/50 ecommended for VarioPrep columns with 50 mm ID or more	1 pz
MNL718256	VP Guard column holder 16 mm (short) VarioPrep guard column holder (short) for 16 mm ID guard column cartridges VP 10/16 recommended for VarioPrep columns with 16 or 21 mm ID	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

