




Chiusure a vite N 9, pronte assemblate

Macherey Nagel

Chiusure a vite N 9, pronte assemblate

Tappo
Setto
Durezza
Spessore mm



Code	Description	Packaging
MNL702080	Chiusure a vite N 9, pronte assemblate, Tappo nero, con foro centrale, Setto Gomma Rossa / FEP incolore, Durezza 45° shore A, Spessore 1 mm	100 pz
MNL702081	Chiusure a vite N 9, pronte assemblate, Tappo rosso, con foro centrale, Setto Gomma Rossa / FEP incolore, Durezza 45° shore A, Spessore 1 mm	100 pz
MNL702082	Chiusure a vite N 9, pronte assemblate, Tappo verde, con foro centrale, Setto Gomma Rossa / FEP incolore, Durezza 45° shore A, Spessore 1 mm	100 pz
MNL702147	Chiusure a vite N 9, pronte assemblateyw, RR or/FEP, 45°, 1.0 N 9 PP screw cap, yellow, center hole Red Rubber/FEP colourless Hardness: 45° shore A Thickness: 1.0 mm	100 pz
MNL702161	Chiusure a vite N 9, pronte assemblate, PP screw cap, blue, center hole	100 pz
MNL702287	Chiusure a vite N 9, pronte assemblate, tr, Sil w/PTFE r, 45°, 1.0 N 9 PP screw cap, transparent, center hole Silicone white/PTFE red Hardness: 45° shore A Thickness: 1.0 mm	100 pz
MNL7022871	Chiusure a vite N 9, pronte assemblate, bl, Sil w/PTFE r, 45°, 1.0 N 9 PP screw cap, blue, center hole Silicone white/PTFE red Hardness: 45° shore A Thickness: 1.0 mm	100 pz
MNL7022881	Chiusure a vite N 9, pronte assemblate,bl, Sil w/PTFE bl, (-) 55°, 1.0 N 9 PP screw cap, blue, center hole Silicone white/PTFE blue, slit Hardness: 55° shore A Thickness: 1.0 mm	100 pz
MNL702402	Chiusure a vite N 9, pronte assemblate, bl, Sil w/Polyimide or, 1.0 N 9 PP screw cap, blue, center hole Silicone white/Polyimide orange Hardness: 50° shore A Thickness: 1.0 mm	100 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

