



**Accessori per CAMERA A VUOTO Chromabond**

**Macherey Nagel**

**Accessori per CAMERA A VUOTO Chromabond**  
Colonne SPE.

Tipo 3 ml, 200 mg



| Code      | Description                                     | Packaging |
|-----------|-------------------------------------------------|-----------|
| MNL730012 | CHROMABOND C18 ec, Volume 3 ml, Capacità 200 mg | 50 pz     |
| MNL730013 | CHROMABOND C18 ec, Volume 3 ml, Capacità 500 mg | 50 pz     |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

