



Colonne per HPLC preparativa POLYGOPREP

Macherey Nagel

Colonne per HPLC preparativa POLYGOPREP



Code	Description	Packaging
MNL7110011000	POLYGOPREP 60-12 pack of 1000 g in plastic container	1 pz
MNL7110021000	POLYGOPREP 100-12 pack of 1000 g in plastic container	1 pz
MNL7110031000	POLYGOPREP 100-20 pack of 1000 g in plastic container	1 pz
MNL7110041000	POLYGOPREP 300-12 pack of 1000 g in plastic container	1 pz
MNL7110051000	POLYGOPREP 1000-30 pack of 1000 g in plastic container	1 pz
MNL7110061000	POLYGOPREP 1000-50 pack of 1000 g in plastic container	1 pz
MNL711007100	POLYGOPREP 60-12 C8 pack of 100 g in glass container	1 pz
MNL7110071000	POLYGOPREP 60-12 C8 pack of 1000 g in plastic container	1 pz
MNL711008100	POLYGOPREP 60-20 C8 pack of 100 g in glass container	1 pz
MNL7110081000	POLYGOPREP 60-20 C8 pack of 1000 g in plastic container	1 pz
MNL711009100	POLYGOPREP 60-12 C18 pack of 100 g in glass container	1 pz
MNL7110091000	POLYGOPREP 60-12 C18 pack of 1000 g in plastic container	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

