



Colonne per HPLC preparativa POLYGOPREP C18

Macherey Nagel

Colonne per HPLC preparativa POLYGOPREP C18



Code	Description	Packaging
MNL711011100	POLYGOPREP 60-80 C18 pack of 100 g in glass container	1 pz
MNL7110111000	POLYGOPREP 60-80 C18 pack of 1000 g in plastic container	1 pz
MNL711018100	POLYGOPREP 100-12 C18 pack of 100 g in glass container	1 pz
MNL7110181000	POLYGOPREP 100-12 C18 pack of 1000 g in plastic container	1 pz
MNL711019100	POLYGOPREP 100-20 C18 pack of 100 g in glass container	1 pz
MNL7110191000	POLYGOPREP 100-20 C18 pack of 1000 g in plastic container	1 pz
MNL711021100	POLYGOPREP 100-50 C18 pack of 100 g in glass container	1 pz
MNL7110211000	POLYGOPREP 100-50 C18 pack of 1000 g in plastic container	1 pz
MNL711024100	POLYGOPREP 300-12 C18 pack of 100 g in glass container	1 pz
MNL7110241000	POLYGOPREP 300-12 C18 pack of 1000 g in plastic container	1 pz
MNL711025100	POLYGOPREP 300-20 C18 pack of 100 g in glass container	1 pz
MNL7110251000	POLYGOPREP 300-20 C18 pack of 1000 g in plastic container	1 pz
MNL711028100	POLYGOPREP 1000-30 C18 pack of 100 g in glass container	1 pz
MNL7110281000	POLYGOPREP 1000-30 C18 pack of 1000 g in plastic container	1 pz
MNL711029100	POLYGOPREP 1000-50 C18 pack of 100 g in glass container	1 pz
MNL7110291000	POLYGOPREP 1000-50 C18 pack of 1000 g in plastic container	1 pz
MNL711031100	POLYGOPREP 60-20 C18 pack of 100 g in glass container	1 pz
MNL7110311000	POLYGOPREP 60-20 C18 pack of 1000 g in plastic container	1 pz
MNL711032100	POLYGOPREP 100-30 C18 pack of 100 g in glass container	1 pz
MNL7110321000	POLYGOPREP 100-30 C18 pack of 1000 g in plastic container	1 pz



Code	Description	Packaging
MNL711470100	POLYGOPREP 60-30 C8 pack of 100 g in glass container	1 pz
MNL7114701000	POLYGOPREP 60-30 C8 pack of 1000 g in plastic container	1 pz
MNL711480100	POLYGOPREP 60-30 C18 pack of 100 g in glass container	1 pz
MNL7114801000	POLYGOPREP 60-30 C18 pack of 1000 g in plastic container	1 pz
MNL711490100	POLYGOPREP 60-50 C8 pack of 100 g in glass container	1 pz
MNL7114901000	POLYGOPREP 60-50 C8 pack of 1000 g in plastic container	1 pz
MNL711500100	POLYGOPREP 60-50 C18 pack of 100 g in glass container	1 pz
MNL7115001000	POLYGOPREP 60-50 C18 pack of 1000 g in plastic container	1 pz
MNL7115401000	POLYGOPREP 100-30 pack of 1000 g in plastic container	1 pz
MNL711590100	POLYGOPREP 60-130 C18 pack of 100 g in glass container	1 pz
MNL7115901000	POLYGOPREP 60-130 C18 pack of 1000 g in plastic container	1 pz
MNL711720100	POLYGOPREP 300-30 C18 pack of 100 g in glass container	1 pz
MNL7117201000	POLYGOPREP 300-30 C18 pack of 1000 g in plastic container	1 pz
MNL711730100	POLYGOPREP 300-50 C18 pack of 100 g in glass container	1 pz
MNL7117301000	POLYGOPREP 300-50 C18 pack of 1000 g in plastic container	1 pz





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

