



Colonne per HPLC preparativa POLYGOPREP C4

Macherey Nagel

Colonne per HPLC preparativa POLYGOPREP C4



Code	Description	Packaging
MNL711022100	POLYGOPREP 300-12 C4 pack of 100 g in glass container	1 pz
MNL7110221000	POLYGOPREP 300-12 C4 pack of 1000 g in plastic container	1 pz
MNL711023100	POLYGOPREP 300-20 C4 pack of 100 g in glass container	1 pz
MNL7110231000	POLYGOPREP 300-20 C4 pack of 1000 g in plastic container	1 pz
MNL711026100	POLYGOPREP 1000-30 C4 pack of 100 g in glass container	1 pz
MNL7110261000	POLYGOPREP 1000-30 C4 pack of 1000 g in plastic container	1 pz
MNL711027100	POLYGOPREP 1000-50 C4 pack of 100 g in glass container	1 pz
MNL7110271000	POLYGOPREP 1000-50 C4 pack of 1000 g in plastic container	1 pz
MNL711690100	POLYGOPREP 300-30 C4 pack of 100 g in glass container	1 pz
MNL7116901000	POLYGOPREP 300-30 C4 pack of 1000 g in plastic container	1 pz
MNL711700100	POLYGOPREP 300-50 C4 pack of 100 g in glass container	1 pz
MNL7117001000	POLYGOPREP 300-50 C4 pack of 1000 g in plastic container	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

