



**Colonne per HPLC preparativa POLYGOPREP NH2**

**Macherey Nagel**

Colonne per HPLC preparativa POLYGOPREP NH2



| Code          | Description                                              | Packaging |
|---------------|----------------------------------------------------------|-----------|
| MNL711012100  | POLYGOPREP 60-12 NH2 pack of 100 g in glass container    | 1 pz      |
| MNL7110121000 | POLYGOPREP 60-12 NH2 pack of 1000 g in plastic container | 1 pz      |
| MNL711013100  | POLYGOPREP 60-20 NH2 pack of 100 g in glass container    | 1 pz      |
| MNL7110131000 | POLYGOPREP 60-20 NH2 pack of 1000 g in plastic container | 1 pz      |
| MNL711014100  | POLYGOPREP 60-30 NH2 pack of 100 g in glass container    | 1 pz      |
| MNL7110141000 | POLYGOPREP 60-30 NH2 pack of 1000 g in plastic container | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

