



Colonne per HPLC preparativa POLYGOPREP Silica SiOH

Macherey Nagel

Colonne per HPLC preparativa POLYGOPREP Silica SiOH



Code	Description	Packaging
MNL7110331000	POLYGOPREP 100-80 pack of 1000 g in plastic container	1 pz
MNL7110341000	POLYGOPREP 100-130 pack of 1000 g in plastic container	1 pz
MNL7110351000	POLYGOPREP 1000-12 pack of 1000 g in plastic container	1 pz
MNL7110361000	POLYGOPREP 1000-20 pack of 1000 g in plastic container	1 pz
MNL7110371000	POLYGOPREP 60-130 pack of 1000 g in plastic container	1 pz
MNL7112401000	POLYGOPREP 60-20 pack of 1000 g in plastic container	1 pz
MNL7112501000	POLYGOPREP 60-30 pack of 1000 g in plastic container	1 pz
MNL7112601000	POLYGOPREP 60-50 pack of 1000 g in plastic container	1 pz
MNL7112701000	POLYGOPREP 60-80 pack of 1000 g in plastic container	1 pz
MNL7115501000	POLYGOPREP 100-50 pack of 1000 g in plastic container	1 pz
MNL7116101000	POLYGOPREP 300-20 pack of 1000 g in plastic container	1 pz
MNL7116201000	POLYGOPREP 300-30 pack of 1000 g in plastic container	1 pz
MNL7116301000	POLYGOPREP 300-50 pack of 1000 g in plastic container	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



