



Colonne per HPLC POLYGOSIL, applicazioni analitiche

Macherey Nagel

Colonne per HPLC POLYGOSIL, applicazioni analitiche



Code	Description	Packaging
MNL71101010	POLYGOSIL 60-5 pack of 10 g in glass container	1 pz
MNL711010100	POLYGOSIL 60-5 pack of 100 g in glass container	1 pz
MNL71102010	POLYGOSIL 60-10 pack of 10 g in glass container	1 pz
MNL711020100	POLYGOSIL 60-10 pack of 100 g in glass container	1 pz
MNL71128010	POLYGOSIL 60-7 pack of 10 g in glass container	1 pz
MNL711280100	POLYGOSIL 60-7 pack of 100 g in glass container	1 pz
MNL71130010	POLYGOSIL 60-5 C8 pack of 10 g in glass container	1 pz
MNL711300100	POLYGOSIL 60-5 C8 pack of 100 g in glass container	1 pz
MNL71131010	POLYGOSIL 60-7 C8 pack of 10 g in glass container	1 pz
MNL711310100	POLYGOSIL 60-7 C8 pack of 100 g in glass container	1 pz
MNL71132010	POLYGOSIL 60-10 C8 pack of 10 g in glass container	1 pz
MNL711320100	POLYGOSIL 60-10 C8 pack of 100 g in glass container	1 pz
MNL71133010	POLYGOSIL 60-5 C18 pack of 10 g in glass container	1 pz
MNL711330100	POLYGOSIL 60-5 C18 pack of 100 g in glass container	1 pz
MNL71134010	POLYGOSIL 60-7 C18 pack of 10 g in glass container	1 pz
MNL711340100	POLYGOSIL 60-7 C18 pack of 100 g in glass container	1 pz
MNL71135010	POLYGOSIL 60-10 C18 pack of 10 g in glass container	1 pz
MNL711350100	POLYGOSIL 60-10 C18 pack of 100 g in glass container	1 pz
MNL71136010	POLYGOSIL 60-5 NH2 pack of 10 g in glass container	1 pz
MNL711360100	POLYGOSIL 60-5 NH2 pack of 100 g in glass container	1 pz



Code	Description	Packaging
MNL71137010	POLYGOSIL 60-10 NH ₂ pack of 10 g in glass container	1 pz
MNL711370100	POLYGOSIL 60-10 NH ₂ pack of 100 g in glass container	1 pz
MNL71138010	POLYGOSIL 60-5 CN pack of 10 g in glass container	1 pz
MNL711380100	POLYGOSIL 60-5 CN pack of 100 g in glass container	1 pz
MNL71139010	POLYGOSIL 60-10 CN pack of 10 g in glass container	1 pz
MNL711390100	POLYGOSIL 60-10 CN pack of 100 g in glass container	1 pz
MNL71142010	POLYGOSIL 60-5 N(CH ₃) ₂ pack of 10 g in glass container	1 pz
MNL711420100	POLYGOSIL 60-5 N(CH ₃) ₂ pack of 100 g in glass container	1 pz
MNL71143010	POLYGOSIL 60-10 N(CH ₃) ₂ pack of 10 g in glass container	1 pz
MNL711430100	POLYGOSIL 60-10 N(CH ₃) ₂ pack of 100 g in glass container	1 pz
MNL71151010	POLYGOSIL 100-5 pack of 10 g in glass container	1 pz
MNL711510100	POLYGOSIL 100-5 pack of 100 g in glass container	1 pz
MNL71152010	POLYGOSIL 100-7 pack of 10 g in glass container	1 pz
MNL711520100	POLYGOSIL 100-7 pack of 100 g in glass container	1 pz
MNL71153010	POLYGOSIL 100-10 pack of 10 g in glass container	1 pz
MNL711530100	POLYGOSIL 100-10 pack of 100 g in glass container	1 pz
MNL71156010	POLYGOSIL 100-5 C ₁₈ pack of 10 g in glass container	1 pz
MNL711560100	POLYGOSIL 100-5 C ₁₈ pack of 100 g in glass container	1 pz
MNL71157010	POLYGOSIL 100-7 C ₁₈ pack of 10 g in glass container	1 pz
MNL711570100	POLYGOSIL 100-7 C ₁₈ pack of 100 g in glass container	1 pz



Code	Description	Packaging
MNL71158010	POLYGOSIL 100-10 C18 pack of 10 g in glass container	1 pz
MNL711580100	POLYGOSIL 100-10 C18 pack of 100 g in glass container	1 pz
MNL71160010	POLYGOSIL 300-7 pack of 10 g in glass container	1 pz
MNL711600100	POLYGOSIL 300-7 pack of 100 g in glass container	1 pz
MNL71168010	POLYGOSIL 300-7 C4 pack of 10 g in glass container	1 pz
MNL711680100	POLYGOSIL 300-7 C4 pack of 100 g in glass container	1 pz
MNL71171010	POLYGOSIL 300-7 C18 pack of 10 g in glass container	1 pz
MNL711710100	POLYGOSIL 300-7 C18 pack of 100 g in glass container	1 pz
MNL71189010	POLYGOSIL 1000-7 pack of 10 g in glass container	1 pz
MNL711890100	POLYGOSIL 1000-7 pack of 100 g in glass container	1 pz
MNL71199110	POLYGOSIL 1000-7 C4 pack of 10 g in glass container	1 pz
MNL711991100	POLYGOSIL 1000-7 C4 pack of 100 g in glass container	1 pz
MNL71199210	POLYGOSIL 1000-7 C18 pack of 10 g in glass container	1 pz
MNL711992100	POLYGOSIL 1000-7 C18 pack of 100 g in glass container	1 pz





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

