



Precolonne fenil-sil, medio-bassa polarità

Macherey Nagel

Precolonne fenil-sil, medio-bassa polarità



Code	Description	Packaging
MNL72310810	Deactivated capillary column without cage Phenyl-Sil deactivated ID: 0.25 mm, length: 10 m	1 pz
MNL72310825	Deactivated capillary column without cage Phenyl-Sil deactivated ID: 0.25 mm, length: 25 m	1 pz
MNL72334810	Deactivated capillary without cage Phenyl-Sil deactivated ID: 0.32 mm, length: 10 m	1 pz
MNL72334825	Deactivated capillary without cage Phenyl-Sil deactivated ID: 0.32 mm, length: 25 m	1 pz
MNL72356010	Deactivated capillary without cage Phenyl-Sil deactivated ID: 0.53 mm, length: 10 m	1 pz
MNL72356025	Deactivated capillary without cage Phenyl-Sil deactivated ID: 0.53 mm, length: 25 m	1 pz

