



Sonda Luminaris color probe High Rox qPCR Master Mix

Thermo Scientific

Luminaris color probe High Rox qPCR Master Mix 250x20ul



Code	Description	Packaging
FL50K0331	Luminaris color probe High Rox qPCR Master Mix250x20ul	250x20ul (*)



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

