



Luminaris color Hi Green qPCR master Mix

Thermo Scientific

Luminaris color Hi Green qPCR master Mix1250x20ul



Code	Description	Packaging
FL50K0393	Luminaris color Hi Green qPCR master Mix1250x20ul	1250x20ul (*)



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

