



Standard per ICP-MS, alluminio, 1000 µg/mL (1000 ppm) in 2-5% acido nitrico (HNO₃)

Reagecon

Standard per ICP-MS, alluminio, 1000 µg/mL (1000 ppm) in 2-5% acido nitrico (HNO₃).



Code	Description	Packaging
LLG04011734	Standard per ICP-MS, alluminio, 1000 µg/mL (1000 ppm) in 2-5% acido nitrico (HNO ₃)	100 ml



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

