



Standard per ICP-MS, alluminio, 1000 µg/mL (1000 ppm) in 2-5% acido nitrico (HNO)

Reagecon

Standard per ICP-MS, alluminio, 1000 µg/mL (1000 ppm) in 2-5% acido nitrico (HNO).



Code	Description	Packaging
LLG04011734	Standard per ICP-MS, alluminio, 1000 µg/mL (1000 ppm) in 2-5% acido nitrico (HNO <sub>3</sub> )	100 ml



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

