



Densimetro ASTM tipo S508HL

Ludwig Schneider GmbH &
Co.KG

Densimetro ASTM S508HL, $-20...+65^{\circ}\text{C}$ in 1°C , 1000-1050 in $0,5 \text{ kg/m}^3$, 380 mm, con certificato DAKKS.



Code	Description	Packaging
LLG04685156	Densimetro ASTM S508HL, '-20...+65°C in 1°C, 1000-1050 in 0,5 kg/m ³ , 380 mm, con certificato DAKKS	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

