



Filtri a membrana isopore 0,22  $\mu\text{m}$ , 25 mm in policarbonato, bianco idrofilo, lisci

Millipore

Filtri a membrana isopore 0,22  $\mu\text{m}$ , 25 mm in policarbonato, bianco idrofilo, lisci.



Code	Description	Packaging
LLG06058008	Filtri a membrana isopore 0,22 µm, 25 mm in policarbonato, bianco idrofilo, lisci, confezione da 100 pz.	100 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

