


**Selettore per diametro tubi****Kartell**

Selettore per tubi in ABS verde. Permette la selezione dei tubi in plastica, vetro o gomma secondo il diametro interno ed esterno.

- Diametri interni disponibili: mm 1 / 2 / 2,5 / 3 / 4 / 4,5 / 5 / 6 / 7 / 8 / 9,5 / 10 / 11 / 12 / 13 / 14 / 15 / 16 / 17 / 18 / 19 / 20 / 21 / 22 / 25.
- Diametri esterni disponibili: mm 1,5 / 2 / 2,5 / 3 / 4 / 5,5 / 6 / 7 / 8 / 9 / 10 / 11 / 12 / 13 / 14 / 15 / 16 / 17 / 18 / 21.



Code	Description	Packaging
LLG06229450	Selettore per tubi, in ABS verde, dim. 90x165 mm	1 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

