



Supporto per capsule Petri Ø35mm

MOTIC

Ø35mm Petri dish holder



| Code        | Description   | Packaging |
|-------------|---|-----------|
| LLG06243521 | Supporto per piastre Petri da Ø 35 mm, per serie AE2000 e AE31E | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

