



**Test per ferro Spectroquant**

**Merck**

Iron Test Spectroquant



| Code        | Description   | Packaging |
|-------------|---|-----------|
| LLG06283155 | Test del ferro Spectroquant® Metodo: fotometrico, kit per test chimici UN 3316, 9, II, (E), MR: 0,005- 5,00mg/l, 250 test | 250 test  |



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

