



Lavavetreria da laboratorio automatica PG 8593 AE-WW-AD, acciaio inox, 835x600x600mm,  
Qmax.=400 L/min

Miele

Lavavetreria da laboratorio con dosaggio liquido e asciugatura aria calda DryPlus.  
Collegabile di serie tramite WiFi o Ethernet.





**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

