



Electronical Digital Thermeter ad 15 lab -50...+200:0,1°C, semiconductor sensor stainless steel,  
35x5mm, 500mm PTFE cable heat resistant up

**Amarell**



Code	Description	Packaging
LLG06490740	Termometro digitale elettronico 15 lab	1 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

