

**Siringhe Combitips advanced Box 2.0, epQuality****Eppendorf**

Siringhe Combitips advanced Box 2.0, epQuality.
Per dispensatori a spostamento positivo.
Adatta anche a liquidi viscosi e volatili.



Code	Description	Packaging
LLG06492975	Siringhe Combitips advanced Box 2.0, epQuality	1 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

