



Olio di calibrazione al silicone, viscosità 500 mPa.s, volume 250 ml, temp. 23 °C

Lamy Rheology

Olio di calibrazione al silicone, viscosità 500 mPa.s, volume 250 ml, temp. 23 °C.



Code	Description	Packaging
LLG06522056	Olio di calibrazione al silicone, viscosità 500 mPa.s, volume 250 ml, temp. 23 °C	1 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

