



**Bottiglia con tappo a vite, 200 ml, in PC-CB**

**Herolab**

Bottiglia con tappo a vite, 200 ml, in PC-CB. Dim. diam. 56 x 112 mm. Per centrifuga. Rpm: 10.000.



Code	Description	Packaging
LLG06522461	Bottiglia con tappo a vite, 200 ml, in PC-CB	6 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

