



Inserto di ricambio in PTFE, per tappo a vite, diam. 13 mm, conf. 20

DWK Life Sciences

Inserto di ricambio in PTFE, per tappo a vite, diam. 13 mm.



| Code        | Description   | Packaging |
|-------------|---|-----------|
| LLG06524332 | Pyrex Screwcap Spare Inserts PTFE, 13 mm pack of 20 | 20 pz     |



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

