



Igrometro Wg. blu, con termometro 0...+35 °C, lungo circa 330 mm, range 1,000-1,200/ div. 0,002 g/ml, temp. di riferimento 20°C

**GECO Gering**

Igrometro Wg. blu, con termometro 0...+35 °C, lungo circa 330 mm, range 1,000-1,200/ div. 0,002 g/ml, temp. di riferimento 20°C.



Code	Description	Packaging
LLG07970837	Igrometro con termometro Wg. blu, range 0...+35 °C, lungo circa 330 mm, 1,000-1,200 in 0,002 g/ml, Tp. 20°C	1 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

