



Termocoppia a isolamento minerale NiCri-Ni tipo K, -200...+800°C, secondo DIN EN 60584, cl.1

**Ludwig Schneider**

Termocoppia a isolamento minerale NiCri-Ni tipo K, -200...+800°C, secondo DIN EN 60584, cl.1.



Code	Description	Packaging
LLG09236709	Termocoppia a isolamento minerale NiCri-Ni tipo K, -200...+800°C, secondo DIN EN 60584, cl.1	1 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

