

**Flash-Chromatography kit I, 400x20****Macherey Nagel**

Flash chromatography kit I consists of: glass column 400/20 1 l eluent reservoir, 100 g silica gel 60, 1 m PTFE tubing, 500 g seasand, glass fiber wadding, silanized



Code	Description	Packaging
MNL727450	Flash-Chromatography kit I, 400x20	1 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

