



Colonna Chromab. SiOH, 3 mL, 500 mg, MPS

Macherey Nagel

CHROMABOND MPS columns SiOH volume: 3 mL, content of sorbent: 500 mg
material: PP with PE filter elements pack of 50 CHROMABOND SPE cartridges
for Gerstel-MPS-autosampler - equipped with cap and cannula



Code	Description	Packaging
MNL730073MPS	Chromab. columns SiOH, 3 mL, 500 mg, MPS	50 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

