



Colonna Chromab. C4, 3 mL, 500 mg

Macherey Nagel

CHROMABOND columns C4 volume: 3 mL, content of sorbent: 500 mg material: PP, with PE filter elements pack of 50



| Code | Description | Packaging |
|-----------|-----------------------------------|-----------|
| MNL730227 | Chromab. columns C4, 3 mL, 500 mg | 50 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

