



Colonna Chromab. PS-MIX, 3 mL, 500 mg

Macherey Nagel

CHROMABOND columns PS-Mix (PS-H+ and PS-OH-) volume: 3 mL, content of sorbent: 500 mg material: PP, with PE filter elements pack of 30



Code	Description	Packaging
MNL730394	Chromab. columns PS-MIX, 3 mL, 500 mg	30 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

