



Colonna Chromab. HR-XAW (45 μ m), 1mL, 30mg

Macherey Nagel

CHROMABOND columns HR-XAW (45 μ m, fine grained) weak mixed-mode polymer-based anion exchanger particle size: 45 μ m volume: 1 mL, content of sorbent: 30 mg material: PP, with PE filter elements pack of 30



Code	Description	Packaging
MNL730728P45	Chromab. columns HR-XAW (45 μ m), 1mL, 30mg	30 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

