



Colonna Chromab. HR-XAW, 3 mL, 60 mg,BIG

Macherey Nagel

CHROMABOND columns HR-XAW weak mixed-mode polymer-based anion exchanger volume: 3 mL, content of sorbent: 60 mg material: PP, with PE filter elements packed in 5 bags with 50 pcs. pack of 250



Code	Description	Packaging
MNL730747250	Chromab. columns HR-XAW, 3 mL, 60 mg,BIG	250 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

