

**Colonna Chromab. HR-XAW, 3 mL, 200 mg****Macherey Nagel**

CHROMABOND columns HR-XAW weak mixed-mode polymer-based anion exchanger volume: 3 mL, content of sorbent: 200 mg material: PP, with PE filter elements pack of 30



Code	Description	Packaging
MNL730748	Chromab. columns HR-XAW, 3 mL, 200 mg	30 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

