



**Colonna Chromab. HR-XAW (45  $\mu$ m),3mL,200mg**

**Macherey Nagel**

CHROMABOND columns HR-XAW (45  $\mu$ m, fine grained) weak mixed-mode polymer-based anion exchanger particle size: 45  $\mu$ m volume: 3 mL, content of sorbent: 200 mg material: PP, with PE filter elements pack of 30



| Code         | Description                               | Packaging |
|--------------|---|-----------|
| MNL730748P45 | Chromab. columns HR-XAW (45 µm),3mL,200mg | 30 pz     |



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

