



Colonna Chromab. HR-XA (45 μ m), 3mL, 60mg

Macherey Nagel

CHROMABOND columns HR-XA (45 μ m, fine grained) strong mixed-mode polymer-based anion exchanger particle size: 45 μ m volume: 3 mL, content of sorbent: 60 mg material: PP, with PE filter elements pack of 30



Code	Description	Packaging
MNL730950P45	Chromab. columns HR-XA (45 μ m), 3mL,60mg	30 pz



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

