



Colonna Chromab. HR-XA, 3 mL, 200 mg

Macherey Nagel

CHROMABOND columns HR-XA strong mixed-mode polymer-based anion exchanger volume: 3 mL, content of sorbent: 200 mg material: PP, with PE filter elements pack of 30



| Code | Description | Packaging |
|-----------|--------------------------------------|-----------|
| MNL730951 | Chromab. columns HR-XA, 3 mL, 200 mg | 30 pz |



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

