

**Chromab. Flash DL 25, colonna vuota****Macherey Nagel**

CHROMABOND Flash DL 25 Flash Direct Load cartridges - empty columns for solid injection  
volumen: 45 mL (for app. 25 g silica or 15 g kieselguhr) material: PP incl. 3 PE filter elements  
(one PE filter element inserted, two aside) inlet: female Luer Lock outlet: male Luer Lock pack of  
50



Code	Description	Packaging
MNL732982	Chromab. Flash DL 25, empty columns	50 pz



**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

