



**Beuta in vetro boro 3.3, collo stretto**

**BRAND**

Beuta in vetro boro 3.3, collo stretto - secondo DIN ISO 1773 - con bordo salvagoccia e graduazione.



| Code        | Description   | Packaging |
|-------------|---|-----------|
| LLG04689213 | beuta in vetro - collo stretto - cap. 100 mL - dim. 64x105 h mm   | 10 pz     |
| LLG04689214 | beuta in vetro - collo stretto - cap. 300 mL - dim. 87x156 h mm   | 10 pz     |
| LLG04689215 | beuta in vetro - collo stretto - cap. 500 mL - dim. 105x180 h mm  | 10 pz     |
| LLG04689216 | beuta in vetro - collo stretto - cap. 1000 mL - dim. 131x220 h mm | 10 pz     |
| LLG04689217 | beuta in vetro - collo stretto - cap. 2000 mL - dim. 166x280 h mm | 10 pz     |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

