



**Bicchiere borosilicato 3.3 DURAN®, forma bassa**

**BRAND**

Bicchiere in boro 3.3 forma bassa - secondo DIN 12331/ISO 3819 - con graduazione e beccuccio.



Code	Description	Packaging
LLG04688952	becher forma bassa, borosilicato 3.3 DURAN® - graduato - 50 mL - diam. 42x60 h mm	10 pz
LLG04688955	becher forma bassa, borosilicato 3.3 DURAN® - graduato - 250 mL - diam. 70x95 h mm	10 pz
LLG04688956	becher forma bassa, borosilicato 3.3 DURAN® - graduato - 2 L - diam. 132x185 h mm	10 pz
LLG04689004	becher forma bassa, borosilicato 3.3 DURAN® - NON graduato - 10 mL - diam. 26x35 h mm	10 pz
LLG04689030	becher forma bassa, borosilicato 3.3 DURAN® - graduato - 400 mL - diam. 80x110 h mm	10 pz



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

