

**Imbuto in vetro interno scanalato****BRAND**

Imbuto in vetro Boro 3.3 - interno scanalato - per filtrazione rapida - pareti inclinate di 60° - bordo superiore arrotondato e levigato - gambo corto.



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG04688809 | imbuto in vetro scanalato - Ø est. 70 mm - gambo Ø est. 8x70 h mm | 10 pz |
| LLG04688810 | imbuto in vetro scanalato - Ø est. 80 mm - gambo Ø est. 10x80 h mm | 10 pz |
| LLG04688811 | imbuto in vetro scanalato - Ø est. 100 mm - gambo Ø est. 10x100 h mm | 10 pz |
| LLG04688812 | imbuto in vetro scanalato - Ø est. 150 mm - gambo Ø est. 16x150 h mm | 10 pz |



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

