

**Pipetta tarata Classe AS a 2 tratti****BRAND**

Pipetta tarata classe AS (alta qualità e dosaggio veloce) - a due tratti - conformità certificata - in vetro AR-GLAS - conforme alla DIN EN ISO 648 - scolamento parziale - anelli di taratura e indicazioni in blu - calibrata per scolamento (Ex) - certificato del lotto incluso.



Code	Description	Packaging
LLG04688046	pipetta tarata classe AS - a 2 tratti - 1 mL	12 pz
LLG04688047	pipetta tarata classe AS - a 2 tratti - 2 mL	12 pz
LLG04688048	pipetta tarata classe AS - a 2 tratti - 5 mL	6 pz
LLG04688049	pipetta tarata classe AS - a 2 tratti - 10 mL	6 pz
LLG04689113	pipetta tarata classe AS - a 2 tratti - 3 mL	6 pz
LLG04689114	pipetta tarata classe AS - a 2 tratti - 15 mL	6 pz
LLG04689115	pipetta tarata classe AS - a 2 tratti - 20 mL	6 pz
LLG04689116	pipetta tarata classe AS - a 2 tratti - 25 mL	6 pz
LLG04689117	pipetta tarata classe AS - a 2 tratti - 50 mL	6 pz
LLG04689121	pipetta tarata classe AS - a 2 tratti - 0,5 mL	12 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

