



Master Mix Stilla compatibile EvaGreen®

Stilla Technologies

Le PCR MIX naica® sono sistemi a 2 componenti pronti all'uso (tampone A e tampone B) ottimizzati per Crystal Digital PCR™ basata su marcatura fluorescente con EvaGreen®.



Code	Description	Packaging
R10050	5X naica™ PCR MIX [Medium]	1,25 ml
R10056	5X naica™ PCR MIX 1.50 ml	
R10106	10X naica™ PCR MIX 0.750 ml	



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

