

# Technical Specification

**PRODUCT : DIETHYLENE GLYCOL RE - Pure****CODE : 346300****MOLECULAR WEIGHT :106.12****CAS N° : 111-46-6****FORMULA : (HOCH<sub>2</sub>CH<sub>2</sub>)<sub>2</sub>O**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Colour	APHA	<= 10
Density d <sub>20</sub> /4	-	1.112 - 1.122
Refractive index at 20°C	-	1.4425 - 1.4525
Boiling point	°C	240 - 252
Water (K.F.)	%	<= 0.3
Assay (GC)	%	>= 99

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in