

Technical Specification

PRODUCT : HEPTANE mixture of isomers RS - PESTIPUR - For pesticide analysis

CODE : 446840

MOLECULAR WEIGHT :100.21

CAS N° : 142-82-5

FORMULA : CH₃(CH₂)₅CH₃

TEST	U.M.	SPECIFICATION
Refractive index at 20°C	-	1.387 - 1.407
Water (K.F.)	ppm	<= 50
Residue on evaporation	ppm	<= 5
Colour	APHA	<= 10
Aromatic compounds	ppm	<= 100
GC-ECD.Individual peak (Lindane)	ng/L	<= 3
Ret.range 1,2,4 trichlorobenzene to decachlorobiphenyle	-	-

For chlorinated compounds analysis; extraction of pesticides.

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in