

# Technical Specification

**PRODUCT : Chloroform RE Pure Stabilized with amylene**

**CODE : 528325**

**MOLECULAR WEIGHT :119,38**

**CAS N° : 67-66-3**

**FORMULA : CHCl<sub>3</sub>**

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Colour	APHA	<=10
Identification (I.R.)	-	Positive
Density at 20°C	-	1.478 - 1.492
Refractive index at 20°C	-	1.4461 - 1.4491
Boiling point	°C	61 - 61.5
Acidity (HCl)	ppm	<=10
Residue on evaporation	ppm	<=20
Water (K.F.)	ppm	<= 300
Assay (GC)	%	>= 99.9
Amylene	ppm	<= 60

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.  
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in