

Technical Specification

PRODUCT : Nitric acid 67.5% (42Be) RE - Pure

CODE : 305500

MOLECULAR WEIGHT :63,01

CAS N° : 7697-37-2

FORMULA : HNO₃

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Density d ₂₀ /4	-	1.395 -1.415
Residue on ignition	ppm	<= 50
Chloride	ppm	<= 30
Sulfate	ppm	<= 50
Heavy metals (as Pb)	ppm	<= 10
Fe	ppm	<= 10
Assay (acidimetric)	%	65.0 - 70.0

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in