

Technical Specification

PRODUCT : HEXANE mixture of isomers RE - Pure

CODE : 339850

MOLECULAR WEIGHT :86,18

CAS N° : 110-54-3

FORMULA : CH₃(CH₂)₄CH₃

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Density at 20°C	-	0.665 - 0.675
Refractive index at 20°C	-	1.375 - 1.383
Boiling point	°C	63 - 70
Water (K.F.)	ppm	<= 200
Aromatic compounds	ppm	<= 30
Residue on evaporation	ppm	<= 50
Total sulphur (S)	ppm	<= 5

It is under the customer responsibility to make sure the purchased product is suitable for his use and/or application.
 CARLO ERBA Reagents would not be guilty of any misuse or mishandling of any of its products, occurring potential damages or user hurts in case of in